

Olive Oil Profiling on FoodScreener 400 and Fourier-80 Benchtop

Pedro Lopes
Lea Heintz


May 31, 2022

Outline

- > Introduction
- > Sample Calibration & Sample preparation
- > Oil-Profiling at the Food Screener (400 MHz) & F80
- > Features
- > Unique selling points
- > Products
- > Targeted Market Segments
- > Configuration
- > Pricing and business models
- > Q&A

Introduction

The process behind food profiling solutions



Edible Oil Profiling - from 400 MHz to 80 MHz

The development story began with 400 MHz, and it continues with 80 MHz

FoodScreener
400 MHz




F80 Benchtop NMR
80 MHz



Why using NMR for food profiling?

Reproducibility

- 50 replicate measurements, including sample preparation




- databases can be built up and combined using raw data generated wherever an Oil-Profiling equipment is available

One quantification reference for all oil samples and oil parameters

- NMR is intrinsically quantitative

Raw data re-analysis

- whenever a so far irrelevant oil component comes into the focus of interest, it is already in the NMR raw data, and can be analyzed retrospectively



Highly efficient holistic method

- the sample preparation for NMR analysis leaves the oils' component "fingerprint" unchanged
- the whole oil profile - from "matrix effects" down to single components - is deciphered in one run

Why using "comparably insensitive" NMR for food profiling?


Detecting tiny variations by data comparison requires **reproducibility**, both of the sample preparation process and the instrumentation! This is prerequisite for database build-up and global spectra comparison, independent from time, persons, instruments, and facilities.

NMR spectroscopy ...

- is **extremely reproducible**

30 replicate NMR measurements, including sample preparation


When superimposed, they seem to be only one single spectrum



- enables to detect what you are looking for (**targeted**), and to find something you did not look for, or even what is unknown (**non-targeted**)

Regression by NMR

Quantification by multivariate regression is applied (instead of direct quantification = integration of a compound's signal). Main Advantage of using linear Regression: statistical models can be applied and concentrations interpolated. This can be done even for NMR peaks that are hidden in the forest of NMR signals and cannot be integrated directly.



| Parameter | Unit | Value | EVOO Reference (IOC) | | |
|-------------------|--------------------|--------|----------------------|--------|------|
| | | | min | max | Flag |
| Free acidity | %w/w as oleic acid | 0.21 | - | 0.80 | ● |
| Peroxide value | mEq O2/kg | 7.4 | - | 20.0 | ● |
| K270 | - | 0.11 | - | 0.22 | ● |
| K232 | - | 1.8 | - | 2.5 | ● |
| Delta K | - | 0.0078 | - | 0.0100 | ● |
| Total polyphenols | mg/kg | 170 | - | - | ○ |
| Linoleic acid | %m/m methyl esters | 7.9 | 2.5 | 21.0 | ● |

The fundament of the quantification of single components or sum parameters by regression is a **database**, containing

- as many as possible **NMR spectra** (both 400 and 80 MHz)
- sample-assigned values of all **parameters** to be quantified, determined by orthogonal (in our case non-NMR) analytical methods (like wet chemistry, HPLC, titration, ...)

Statistical processing of NMR data then generates regression models for each parameter by correlating specific variations with metadata values.

Then, these regression models are used to quantify corresponding parameters in samples to be analyzed.

Introduction

Main instrumental features

- method installation via an [Installer](#) (analogue to all other FoodScreening methods)
- experiments/reports fully automated for both platforms (SampleTrack)
 - sample setup via [Oil Easy Dialog](#) menu
 - F80 - sample-ID by SampleTrack
 - FoodScreener - barcode labels
 - automated [Analysis Report](#) generation
- F80/FoodScreener usable with or without Sample Changer
- [Standardized Work Instructions](#) for oil sample preparation available
- Oil-Profiling 1.0 Solution [Manuals](#) will be available soon



The Oil-Profiling experimental set

80 MHz platform

- standard ^1H -NMR spectrum, with parameters optimized for oil samples ("zg" pulse sequence)

400 MHz platform

- standard ("zg") ^1H -NMR spectrum, with parameters optimized for oil samples
- ^1H -NMR spectrum where intense lipid signals are suppressed - a "noesy" pulse sequence containing a shaped pulse usable for signal suppression is used
- ^1H -JRES experiment for further spectral information (e.g. easier access to signal multiplicity in regions with overlapping signals)

total experimental time

12 minutes



25 minutes (including sample temperature stabilization, atma, lock, shim, and pulsecal)

Pipette Calibration

Sample Preparation

Pipette Calibration


- The use of positive displacement pipettes (PosD) is necessary because of the low viscosity of the solvent CDCl_3
- Pipettes are initially calibrated according to Bruker Work Instructions. Calibration is done by identifying the correct Oil/ CDCl_3 ratio. This is done automatically by NMR. Typically the volume of oil to be pipetted will be $\sim 214 \mu\text{l}$
- The use of USC caps for the NMR tubes (the purple ones) is necessary to prevent solvent evaporation

Positive Displacement Pipette



olive oil samples in CDCl_3 situation 14 day after sample preparation

Pipette Calibration Instructions




Sample Preparation & Quality Control Measures

- Sample preparation is the same for the F80 and the 400 MHz Food Screener & according to the Bruker Instructions (Video at the end of the presentation)
- Quality Control is done with reference samples from Bruker (not Oil) and with one Oil sample from the customer (can be any oil sample chosen by the customer).
- 400 MHz: 3 reference samples required + 1 Oil reference sample
- 80 MHz: 1 reference samples required + 1 Oil reference sample

| AV400 Reference samples | Produced/prepared by | Further informations or explanations |
|---|--|---|
| QuantRef | Bruker BioSpin GmbH (Germany) | validated solution for quantification |
| Sucrose | Bruker BioSpin AG (Switzerland) | certified* product "2mM sucrose in the NMR tube" |
| TempCal | Bruker BioSpin AG (Switzerland) | certified* product "99.8% MeOD in the NMR tube" |
| F80 Benchtop Reference samples | Produced/prepared by | Further informations or explanations |
| QuantRefB | User (prepare according to flow-chart in Chapter 9.3) | validated solution for quality control and quantification |
| Instrument-independent Reference sample | Produced/prepared by | Further informations or explanations |
| OilRef | User (prepare according to flow-chart in Chapter 7, as any other oil sample) | <ul style="list-style-type: none"> define ~200 ml of any olive oil as "oil for Oil-Ref", and mark it accordingly when consumed, continue likewise with a further 200 ml portion of olive oil prepare the OilRef sample in the course of the preparation of other oil samples one OilRef samples for all NMR instruments |

* Reference sample has to be ordered *always together* with the associated certificate

Sample Preparation Instructions




"Evolution" of the oil sample preparation method

original SOP (weight-based, manual)

EXCEL-guided weighing of oil and CDCl_3 amount

keep this ratio as constant as possible

| Oil weight (mg) | Final weight (mg) | weight Oil / weight CDCl_3 | CDCl_3 volume difference to ref (μl) |
|-----------------|-------------------|-------------------------------------|--|
| 140,00 | 1176,00 | 13,51351351 | 0,00 |
| 139,60 | 1172,10 | 13,52058111 | -0,36 |
| 140,90 | 1184,70 | 13,49875455 | 0,77 |



Quantos



with more and more samples, manual weighing becomes a hassle ...

"interim" SOP (weight-based, semi-automated)

final SOP (volume-based, manual)

"F80-compatible" low-cost solution




Oil-Profiling at the Food Screener (400 MHz) & F80

Olive oil spectra - the lipid profile

Some signals always are quasi-identical, (even for oil types other than olive oil)

Some other signals vary within a certain range from olive oil sample to sample



Identical Reports for the F80 & the 400 MHz

Model:

- oil type (olive oil vs. others)
- country of origin (Spain, Italy, Greece)
- Extra virgin vs. other qualities

Analysis Report Olive Oil Profiling 400 MHz™

Sample ID: 108928_210212

Information/Declaration provided by customer:

Type of Oil: Olive Oil
 Country of Origin: Greece
 Quality Grade: undefined
 Harvest Period: undefined

Disclaimer: this information will affect the applicability and validity of analyses and results.

Measuring Date: 16-Feb-2021 16:00:07
 Reporting Date: 02-Mar-2022 16:39:58, 5 pages, Version 1.0.0

Results Summary

| Type of Analysis | Result | Status |
|--------------------------------------|--------------------------------|--------|
| Analysis of Origin | | |
| Country Greece | Consistent | ● |
| Quantitative Analysis | Consistent with IOC references | ● |
| Comparison to Reference Group | Typical | ● |

The data analysis is performed at Bruker BioSpin GmbH (Ettlingen, Germany) according to testing method Olive Oil Profiling 400 MHz. All results solely refer to the tested sample as provided by the customer.

Quantitative Analysis

(Analysis-ID: OO1400-REG-20220302)

| Parameter | Unit | Value | EVOO Reference (IOC) | | |
|-----------------------|------------------------|--------|----------------------|--------|------|
| | | | min | max | Flag |
| Free acidity | %w/w as oleic acid | 0.21 | - | 0.80 | ● |
| Peroxide value | mEq O2/kg | 7.4 | - | 20.0 | ● |
| K270 | - | 0.11 | - | 0.22 | ● |
| K232 | - | 1.8 | - | 2.5 | ● |
| Delta K | - | 0.0078 | - | 0.0100 | ● |
| Total polyphenols | mg/kg | 170 | - | - | ○ |
| Linoleic acid | %m/m methyl esters | 7.9 | 2.5 | 21.0 | ● |
| Linolenic acid | %m/m methyl esters | 0.80 | - | 1.00 | ● |
| Oleic acid | %m/m methyl esters | 75.4 | 55.0 | 83.0 | ● |
| Palmitic acid | %m/m methyl esters | 11.2 | 7.5 | 20.0 | ● |
| Palmitoleic acid | %m/m methyl esters | 1.0 | 0.3 | 3.5 | ● |
| Stearic acid | %m/m methyl esters | 2.2 | 0.5 | 5.0 | ● |
| Wax content | mg/kg | 7 | - | 150 | ● |
| Erythrodiol + Uvaol | % total sterols | 1.9 | - | 4.5 | ● |
| b-Sitosterol apparent | % total sterols | 94.7 | 93.0 | - | ● |
| Total sterols | mg/kg | 1321 | 1000 | - | ● |
| Total MUFA | % of total fatty acids | 75.9 | - | - | ○ |
| Total PUFA | % of total fatty acids | 8.8 | - | - | ○ |
| Total TFA | % of total fatty acids | 0.07 | - | - | ○ |
| Total SFA | % of total fatty acids | 13.7 | - | - | ○ |

MUFA mono-unsaturated fatty acids
 PUFA poly-unsaturated fatty acids
 TFA trans unsaturated fatty acids
 SFA saturated fatty acids

Video: Oil Profiling F80 & 400 MHz



Please Watch the video starting at minute 26 from the full recorded session: [LINK](#)

Features



Features

| | Olive Oil Profiling 1.0 | Olive Oil Profiling 2.0 |
|---|-------------------------|-------------------------|
| Compliance check for declared country of origin | Spain, Greece | Addition of Italy |
| Compliance check for Extra Virgin OO | Not available | Yes |
| Quantitative Analysis of IOC regulated parameters | 15 parameters | 15 parameters |
| Detection of atypical profiles | Yes | Yes |

- Fully automated analysis
- Rapid analysis: 12mn/sample on Fourier 80 vs 25mn/sample on FoodScreener
- Fast sample preparation: roughly 3mn/sample.
- Under work: ISO17025 accreditation

Unique Selling Points



Unique selling points

Overview of competing techniques

| | Block Chain Genetics (DNA) | Panel Test (Sensory) | Wet Chemistry | GC or LC FID/Fluores/MS | UV-VIS NIR | NMR |
|-------------------------------|----------------------------|-------------------------------|----------------------|-----------------------------|----------------|--------------------------|
| Authenticity | Country of Origin | | | | | Country of Origin |
| | Cultivar | | | | | Cultivar |
| | | | | Blends/dilutions | | Blends/dilutions |
| Quality Grade | | EVOO/VOO Attributes & Defects | | | | EVOO/VOO Screening |
| Quality Phys-Chemistry | | | Tritation Gravimetry | GC-FID for Fatty Acids E.E. | UV 270/232 NIR | Most relevant parameters |

- Combine analysis of country of origin, quality grade and chemical parameters, with one single measurement !
- Screening for quality grade with NMR reduces the number of samples to be tested by panel test (sensory), as this requires experts, and can not exceed 10 tests/day to be reliable.
- Fully automated and easy to use methods
- Fast analysis
- No specialist required to operate the system.
- Ready to use method: no method development required !

Products

Product families



Fourier 80 with automation



FoodScreener

Targeted Market Segments and Value Proposition



Targeted Market Segments

Preferred solution*

| | | Fourier 80 | FoodScreener |
|-------------------------------------|-------------------------------------|------------|--------------|
| Olive Oil Bottlers | | | |
| Commercial Service Providers | Main Hub | | |
| | (International) Subsidiaries | | |
| | Laboratories specialized only in OO | | |
| Governmental laboratories | | | |

* Preferred / recommended product in green. Nevertheless, upon need and request of customer, both products can be sold without restriction to all market segments.

Value Proposition: Olive Oil Bottlers & Valorizers



- = A compact, benchtop solution, easy to implement and to operate !
- = Saving time and money in insourcing analysis (Fast ROI)
- = Replace other techniques like e.g. NIR
- = Quality & authenticity control of olive oil at both purchase and selling side
- = **Brand protection**
- = Use NMR as **marketing tool** to sell premium Olive Oils with tested origin and quality

Value Proposition: Commercial Testing Labs



- = Expand product portfolio & sales revenues with innovative new method
- = Differentiate from competitors
- = Big labs can equip subsidiaries with a benchtop Fourier 80 system
- = Fast ROI, low TCO



Value Proposition: Public & Government Accounts Official Control, Customs & Border Protection, Regulatory Bodies



- = Check if olive oils are compliant with label information & regulation
- = Fight fraud on origin, ensure fair competition
- = New business model (Flat Rate) offer more flexibility and allow to measure more samples
- = Possibilities of cooperation to further extend the method
- = Possibility to access method details for validation purposes or in case of prosecution.

Configuration



Configuration

| | FoodScreener | | Fourier 80 | |
|------------------|---|---------------------------------|--|------------------------------|
| | Options | Price (EUR)* | Options | Price (EUR)* |
| Platform | AV4400FOOD 1H inverse probehead, Boss3 shimsystem, BCU, TopSpin and SampleTrack.. MSASC400SBAIC magnet AH0070 transfer line | 456,898 151,768 3,755 | UH0085 1H/13C system Sample temperature of 25°C TopSpin and SampleTrack Bench PC, touchscreen, cleaning kit, etc. Installation Support included (2 hours of phone) | 83,104 |
| Sample changer | Included in platform 60 pos | | FH0040 Optional (60 pos) | 30,962 (ACA) 32,954 (IND) |
| Barcode reading | Automatical Samples can be put on any place on the sample changer | | Manual Samples need to be placed on a specific position on the sample changer | |
| Olive Oil Module | AH0604 | 13,433 | UH0085OL | 10,000 |
| | Remark: no pH adj. | | | |

* Prices are indicative and may change from time to time, as well as region-dependent. Current prices can be found in SAP

Pricing and business models



Olive Oil Profiling Analyses Options

| | FoodScreener | | Fourier 80 | |
|------------------|---|----------|---|---------------|
| | Options | Segment | Options | Segment |
| Pay-per-use | Not available | | Not available | |
| Pre-paid bundles | 500 bundle LYOO_0500 (13000€) | ALL | 500 bundle, 12 months LYOO_0500_80 (10000€) | ALL |
| | 1000 bundle LYOO_1000 (22000€) | ALL | 1000 bundle, 12 months LYOO_1000_80 (16000€) | ALL |
| Annual Flatrates | Capped to 10000, all food LY_FR1Y_GOV (45000€) | GOV only | Capped to 1000/y LYOO_FR1Y_GOV (8000€) | GOV only |
| | Capped to 600, all food LY_600_1Y_GOV (20000€) | GOV only | Capped to 1000/y LYOO_FR1Y_BOT (8000€) | Bottlers only |
| | | | Quant only LYOO_QFR1Y_BOT (4000€) | Bottlers only |

* Prices are indicative and may change from time to time, as well as region-dependent. Current prices can be found in SAP

Olive Oil Profiling Analyses Options

- **Adjusted Business Model for Gov Accounts (FoodScreener)**
 - Gov accounts = official control laboratories, customs laboratories and regulatory bodies only. **Not valid for universities** and other ACA customers, even if the lab is owned by government.
 - The annual Flat Rate Model is valid for any food matrix
 - The annual Flat Rate models, offer more flexibility to Government accounts and incentive to test more (the more samples analyzed, the cheaper the cost per sample)
 - The flatrates are valid **for 1 calendar year**
 - With the purchase of the first flatrate, the customer can order AH0610_GOV, which offers the following items complimentary:
 - Installation of parameter sets for any missing food module(s); e.g. Olive Oil
 - Free analyses till the end of the calendar year (limited to 600 max)

Q&A

Q&A

Q: How is the Quantification done (e.g. linear regression)?

A: We do not quantify any oil component directly at the moment. All the components are quantified using linear regression. In other words, for every component we quantify, we have a collection of corresponding concentration ranges obtained by non-NMR methods (e.g. HPLC, wet chemical analysis, etc.) to which we apply a linear regression method in order to get the concentration of a given oil component for a given oil sample.

Q: How is peak deconvolution done? A: It is not done so far.

Q: How are the limit ranges (LOQ) obtained?

A: LOQ was not determined so far because the limit of detection (LOD) was not determined so far. Typically, $LOQ = LOD \times 2/3\sigma$. In order to obtain LOD a dilution series is performed. Spiking experiments can also be performed.

Q: Can Total Cholesterol quantification be performed?

A: not automatically with the current version of the profiler. However, if the signal(s) are sufficiently isolated, quantification could be easily done by performing manual integration and using the Eretic method. If the signal is overlapping with other signals, quantification could still be possible to some degree by performing signal deconvolution by peak fitting.

Q: Can Total Fatty Acids and not only Olive Oils be quantified by this method?

A: This is possible but not without adjustments to the current method: so far we quantify by linear regression and this method always requires a database in the background. Presently we have only Olive Oils in the database and regression models for Greek, Spanish & Italian. Direct quantification of any component always requires a database. Direct quantification is a difficult process to establish for fatty acids. Which often have overlapping signals at 400 MHz. Another way is to perform lineshape analysis but this was also not implemented so far.


Q&A


Q: What is the difference between the F80 and the 400 MHz? Will the F80 deliver less accurate results?

The first difference between both is the sensitivity: if all the conditions are the same (same sample concentration, quality, etc) you gain a factor of 11 in S/N ratio by going from the 80 MHz to the 400 MHz. This is established according to the equation:

$$\frac{S/N_{80}}{S/N_{400}} = \left(\frac{400}{80}\right)^{3/2} = 11,2$$

The second difference is that the resolution is reduced: signals that are well separated at 400 MHz may be partially or fully overlapping at 80 MHz. Nevertheless, most informations are available – this is not obvious when visually inspected, but "hidden" e.g. in lineshapes ... Our statistical methods can pick out them also from an 80 MHz spectrum.





Innovation with Integrity